

Beamline Macros

Change Beamline Apertures and Move Diagnostic Stages

- **setBL** *beamline_component value*

Options for *beamline_component*:

1. **vahgap** - variable aperture horizontal gap
 - a. *value* = gap in mm
2. **vavgap** - variable aperture vertical gap
 - a. *value* = gap in mm
3. **esgap** - exit slit gap
 - a. *value* = gap in um
4. **bds3** - beam diagnostic stage 3 element
 - a. *value* = out, yag, OR au
5. **bds4** - beam diagnostic stage 4 element
 - a. *value* = out, yag, au, OR pd

- **getBL** *beamline_component*

Options for *beamline_component*:

1. **vahgap** - variable aperture horizontal gap
2. **vavgap** - variable aperture vertical gap
3. **esgap** - exit slit gap
4. **bds3** - beam diagnostic stage 3 element
5. **bds4** - beam diagnostic stage 4 element

Open and Close Photon Shutters and Gate Valves

- **setPSH** *state* - open or close beamline shutter

Options for state:

1. open OR close

- **statPSH** - Displays state of photon shutter.
- **setVVR** *state* - open or close endstation gate valve

Options for state:

1. open OR close

- **statVVR** - Displays state of photon shutter.

Turn On and Off Silicon Drift Detector

- **setSDD** *state*

Options for state:

1. on OR off

Monochromator Macros

- **getE** - Display current energy of the monochromator
- **setE** *energy*
energy in eV at the current monochromator position

- **moveE** *energy*

Options for *energy*:

1. Ni LEG: 150 - 900 eV
2. Au LEG: 95 - 900 eV
3. Au HEG: 300 - 3000 eV

Notes: Moves energy to precisely the specified value. When moving down the energy overshoots and comes back. For **EScan**, a **moveE** is not necessary as this is executed in all cases.

- **moveG** *grating*

Options for *beamline_component*:

1. NiLEG
2. AuLEG
3. AuHEG

- **moveM** *mirror*

Options for *beamline_component*:

1. Nickel
2. Carbon
3. Silicon
4. Gold

Undulator Macros

- **setEPU** *epu_component value*

Options for *beamline_component*:

1. **polar** - polarization

Options for *value*:

- a. cl - circular left
- b. cr - circular right
- c. lh - linear horizontal
- d. lvp - linear vertical positive
- e. lvn - linear vertical negative
- f. li - linear inclined (include angle after)

2. **harmonic** - harmonic

- a. *value* = 1, 2, 3, 4, OR 5

3. **offset** - EPU gap offset

- a. *value* = offset in mm

Current Amplifier Macros

- **statAMP** - Displays the current sensitivity of the current amplifiers
- **setAMP** *amp sens*

Options for *amp*:

1. *a* - I0 amplifier
2. *b* - TEY amplifier

Options for *sens*:

- a. *up* OR *down* - relative change
- b. *1pA/V* OR *2nA/V* - absolute

Notes: The exhaustive list of gains are 1, 2 and 5 with a range of 1 mA/V to 1 pA/V and all orders of magnitude in between.

Beam Status Macros

- **chkBEAM** *parameter state*

Options for *parameter, state* is on OR off:

3. ring - checks ring state and safety shutters
4. psh - checks state of shutters
5. vvr - checks state of gate valves
6. on OR off - turns all on OR off
7. status - display status of checking

- **setBEAM** *state* - ensures beam is available

Options for *state*:

1. on OR off

Special Scanning Macros

- **Escan** *start end1 intervals end2 intervals2 ... time [fixQ]*

Notes: *fixQ* is only available for four.

- **Tscan** *time_interval total_time*

Standard Motor Macros

- **ascan** *mot start end intervals time*
- **rscan** *mot start end1 intervals end2 intervals2 ... time*
- **lup** *mot rel_start rel_end intervals time*
- **mesh** *mot1 start end intervals mot2 start end intervals time*
- **umv** *mot1 position1 mot2 position etc*
- **umvr** *mot1 rel_position1 mot2 rel_position etc*

RIXS Macros

Sample Stage Macros

- **setSS** *user preset_num desc* - sets to current position
Options for *preset_num*:
 1. *1* through *5*Options for *desc*:
 1. A single string, needs to be connected
- **setSS** *user preset_num desc abs_ssh abs_ssd abs_ssv abs_ssa*
- sets to the position to current position
Options for *preset_num*:
 1. *1* through *5*Options for *desc*:
 1. A single string, needs to be connected*abs_ssh*, *abs_ssd*, *abs_ssv* and *abs_ssa* are the positions of the sample stage motors.
- **moveSS** *preset OR user user_preset* - moves to preset position
Options for *preset*:
 1. *tran* - transfer position
 2. *meas* - measure position
 3. *diode* - photodiode
 4. *cryo* - screw for cryoThe flag *user* is used explicitly and the numbers refers to user preset that was set in **setSS** and can be checked with **statSS**
- **statSS** - checks the current absolute position of the sample stage and displays the user defined presets.

Spectrometer Macros

- **setXES** *option value*

Options for *option*:

1. **grating** - grating for XES spectrometer

Options for *value*:

1. **XLEG**
 2. **LEG**
 3. **MEG**
 4. **HEG**
2. **pinhole** - pinhole for XES spectrometer (no value)
 3. **mask** - mask for XES spectrometer (no value)
- **moveXES** *energy tilt_offset*
Options for *energy*:
 1. Acceptation energy in eV with grating rangeOptions for *tilt_offset*:
 1. Tilt offset in degrees, usually -2 to 4 degrees